

Silicone Simulation Pad Model for Medical Teaching and Practice HS Code: 9023009000

1. Usage: this silicone simulation pad models help to prep up for a real-life scenario, designed to practice and enhance the skill sets of medical, dental, veterinary students.
2. Type: suture pad, IV infusion pad, injection pad, laparoscopic suture pad, dental pad, acupuncture pad.
3. Material: made of high quality silicone, elastic and smooth, ensure that its texture is close to real skin with a standard approved thickness.
4. Durable material, fine workmanship. Can reuse after gently remove the suture.
5. Easy to carry. Make the teaching and practice convenient.

Suture Pad Model

Usage: practice for wound cut, suture, knot, remove suture.

Features: clear wound, simulate real skin structure of three layer –skin layer, fat layer and muscle layer. Inner protective mesh for choose.

Product name: silicone simulation wound suture pad model

Reference no: FH2311001

Product size: 170*120*12 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311002

Product size: 170*120*12 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311003

Product size: 170*120*12 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Product name: silicone simulation wound suture pad model


Reference no: FH2311004

Product size: 200*150*13mm

Package: nonwoven bag individually,
one piece in a box, 40 pcs/ctn

Carton size: 47.5*33.5*31cm

G.W.: 15.5 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311005

Product size: 200*200*13 mm

Package: nonwoven bag individually,
one piece in a box, 30 pcs/ctn

Carton size: 48.5*47.5*24 cm

G.W.: 16 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311006

Product size : 155*115*20 mm

Package: PE bag individually,
one piece in a box, 60 pcs/ctn

Carton size: 41.5*34.5*20 cm

G.W.: 15 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311007

Product size : 155*115*20 mm

Package: PE bag individually,
one piece in a box, 60 pcs/ctn

Carton size: 41.5*34.5*20 cm

G.W.: 14 kgs



Product name: silicone DIY wound suture pad model

Reference no: FH2311008

Product size with box: 175*110*27 mm(box)

Package: one piece in a box, 60 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 15 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311009

Product size with box: 170*120*12 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Product name: silicone DIY wound suture pad model

Reference no: FH2311010

Product size with box: 150*110*13 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 11 kgs



Product name: silicone simulation wound suture pad model

Reference no: FH2311011

Product size with box: 150*110*13 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 11 kgs



3D Suture Pad Model

Usage: design for laparoscopic soft tissue suture, to practice wound cut, suture, knot, remove suture.

Features: have different types of wound for practice. Simulate real skin structure to improve the surgical skill.

Product name: silicone 3D simulation wound suture pad model

Reference no: 3DFH2311001

Product size: 150*150*15 mm

Package: nonwoven bag individually,
one piece in a box, 60 pcs/ctn

Carton size: 53*34*33 cm

G.W.: 11 kgs



Product name: silicone 3D simulation wound suture pad model

Reference no: 3DFH2311002

Product size: 140*140*15 mm

Package: nonwoven bag individually,
one piece in a box, 60 pcs/ctn

Carton size: 53*34*33 cm

G.W.: 11 kgs



Dental Suture Pad Model

Usage: design for gingival wound cut, suture, knot, remove suture.

Features: simulate different types of gingival wound for different types of suture, such as continuous suture, intermittent suture, lockstitch suture, horizontal mattress suture, corner suture, 8-shaped intermittent suture, suspension suture, anchor suture.

Product name: silicone 3D simulation wound suture pad model

Reference no: YYFH2311001

Product size: 170*122*18 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 10 kgs



Product name: silicone 3D simulation wound suture pad model

Reference no: YYFH2311002

Product size: 170*122*18 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 12 kgs



IV infusion Pad Model

Usage: design for practice injection and infusion.

Features: simulate real skin structure of three layer –skin layer, fat layer, muscle layer and including 4 sealed blood vessel. After injection, there is blood return. For repeat use.

Product name: silicone IV infusion pad model


Reference no: ZS1203001

Product size: 170*100*17 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24cm

G.W.: 13 kgs



Peel and Cut Pad Model

Usage: laparoscopic surgery practice.

Features: simulate real abdominal soft and features to practice together with the laparoscopic.

Product name: silicone peel and cut pad model

Reference no: BL23120001

Product size: 138*138*15 mm

Package: nonwoven bag individually,
one piece in a box, 60 pcs/ctn

Carton size: 53*34*33 cm

G.W.: 11 kgs



Blood Collection Pad Model

Usage: design for practice blood collection and infusion.

Features: with a sealed blood vessel. Before use, evacuate air by a syringe in one end and infuse simulate blood in the other end. There will have blood collection when inject. After use, push the simulate blood back to blood vessel.

Product name: silicone blood collection pad model

Reference no: HX1122001

Product size: 180*110*13 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Product name: silicone blood collection pad model

Reference no: HX1122002

Product size: 115*75*20 mm

Package: nonwoven bag individually,
one piece in a box, 50 pcs/ctn

Carton size: 41.5*34.5*24 cm

G.W.: 13 kgs



Suture practice kit



Product name: large suture kit
Reference no: FHT021
Component: suture pad x1;
 leather case x1;
 mosquito forceps x1;
 Iris scissors x1;
 hegar needle holder x1;
 scalpel handle x1;
 adson toothed forceps x1;
 suture thread x20;
 surgical blades x6;
 latex gloves x2
Package: one set per box,
 40 box per carton
Carton size: 485 x465 x345 mm
Gross weight: 21kgs



Product name: small suture kit
Reference no: FHT031
Component: suture pad x1;
 leather case x1;
 mosquito forceps x1;
 Iris scissors x1;
 hegar needle holder x1;
 scalpel handle x1;
 adson toothed forceps x1;
 suture thread x20;
 surgical blades x6;
 latex gloves x2
Package: one set per box,
 40 box per carton
Carton size: 430 x310 x420 mm
Gross weight: 19kgs



Product name: small suture kit
Reference no: FHT003
Component: suture pad x1;
 leather case x1;
 mosquito forceps x1;
 Iris scissors x1;
 hegar needle holder x1;
 scalpel handle x1;
 adson toothed forceps x1;
 suture thread x5;
 surgical blades x5;
 latex gloves x2
Package: one set per box,
 40 box per carton
Carton size: 430 x310 x420 mm
Gross weight: 18kgs



Product name: large suture kit
Reference no: FHT016
Component: suture pad x1;
 leather case x1;
 mosquito forceps x1;
 Iris scissors x1;
 hegar needle holder x1;
 scalpel handle x1;
 adson toothed forceps x1;
 suture thread x5;
 surgical blades x5;
 latex gloves x2
Package: one set per box,
 40 box per carton
Carton size: 485 x465 x345 mm
Gross weight: 18kgs