

# MALE EXTERNAL CATHETERS

## **NANCHANG YILI MEDICAL INSTRUMENT CO., LTD**

Office: 101# Building, Greenland International Exp Center, Jiulong Road 1177, Honggutan District, Nanchang, Jiangxi, China

Factory: Jinxian Medical Industrial Park, Nanchang, Jiangxi, China

Tel: +86-791-88252326 Fax: +86-791-88251326 Post Code: 330046


E-mail: [sales@yilimed.com](mailto:sales@yilimed.com)

[www.yili-medical.com](http://www.yili-medical.com)



[www.yili-medical.com](http://www.yili-medical.com)

# EasyThru®






# Male External Catheters

EasyThru<sup>®</sup>

Male external catheter is a body worn urine drainage devices which may be attached to a vacuum device or a urine collection bag. It is rolled onto the penis and secured by an adhesive strap or self adhesive.



## LATEX MALE EXTERNAL CATHETER

### Basic details:

- 1. Made of medical grade latex.
- 2. Composition: latex external catheter, adhesive strap
- 3. Used for male patients who suffer paralysis or uroclepsia.
- 4. Sterilized by ethylene oxide gas.
- 5. Certificate: ISO13485, CE, FDA
- 6. Shelf time: 3 years

### Advantages:

- 1. External catheters are the least invasive for draining urine and are available as penile sheath catheters for men
- 2. Reduce the risks of damaging the urethra and infection.



Code	Size
YL-LEC-20	20mm
YL-LEC-25	25mm
YL-LEC-30	30mm
YL-LEC-35	35mm
YL-LEC-40	40mm

# Male External Catheters

## SILICONE MALE EXTERNAL CATHETER

EasyThru®

### Basic details:

1. Made of medical grade silicone;
2. Composition: silicone external catheter, adhesive strap;
3. Sterilized by ethylene oxide gas;
4. Certificate: ISO13485, CE, FDA;
5. Shelf time: 3 years.

### Adhesive strap

#### Advantages:

1. Latex free, all silicone construction provides excellent breath-ability.
2. Available size selection helps ensure a secure fit.
3. Transparent material allows for skin health monitoring while the catheter is in place
4. The catheter provides great safety for both users and caregivers, as it prevents skin

Code	Size	Type	Code	Size	Type
YL-NSEC-25	25mm	Adhesive strap	YL-NSEC-30	30mm	Adhesive strap
YL-NSEC-35	35mm	Adhesive strap	YL-NSEC-40	40mm	Adhesive strap

### Self-adhesive

#### Basic details:


1. Made of medical grade silicone;
2. Sterilized by ethylene oxide gas;
3. Certificate: ISO13485, CE, FDA;
4. Shelf time: 3 years.

#### Advantages:


1. Latex free, all silicone construction provides excellent breath-ability.
2. Available size selection helps ensure a secure fit.
3. Transparent material allows for skin health monitoring while the catheter is in place.
4. Safe, strongly skin friendly adhesive bonding has good initial cohesiveness, and good holding power time.
5. The catheter provides great convenience and safety for both users and caregivers, as it prevents skin.

#### Types:


1. pop on: shorter sheath length;



2. standard: regular length and adhesive length;



3. wide band: full length adhesive for adding security.



Code	Size	Type	Code	Size	Type	Code	Size	Type
YL-SECP-25	25	Pop on	YL-SEC-25	25	Standard	YL-SECW-25	25	Wide band
YL-SECP-29	29	Pop on	YL-SEC-29	29	Standard	YL-SECW-29	29	Wide band
YL-SECP-32	32	Pop on	YL-SEC-32	32	Standard	YL-SECW-32	32	Wide band
YL-SECP-36	36	Pop on	YL-SEC-36	36	Standard	YL-SECW-36	36	Wide band
YL-SECP-40	40	Pop on	YL-SEC-40	40	Standard	YL-SECW-40	40	Wide band



# Male External Catheters

## - Packaging

EasyThru®



( Blister individual )



( Hard blister individual )



(box)



CTN: 50pcs/box, 1000pcs/ctn;  
MEAS: 55\*28\*57CM;  
295000pcs/20GP; 636000pcs/40GP.